

elcometer® 134S Chlor Test for Blast Cleaned Surfaces

Can be used in accordance with: ISO 8502-6, IO DIS 8502-9



Chloride Salts left on the surface before the first coat is applied can result in the coating system being forced off the surface by corrosion or blistering before the full life of the coating has been reached.

To ensure that the chloride has been removed it is essential that the surface is tested before the coating is applied.

Features

- Simple low cost test for field testing of Chloride Ion on surfaces.
- No needles or conductivity meter required.

Specifications

	Elcometer 134A	Elcometer 134S	Elcometer 134W
Measuring Range	1– 50ppm ($\mu\text{g}/\text{cm}^2$)	1 – 60ppm ($\mu\text{g}/\text{cm}^2$)	10 – 2000ppm ($\mu\text{g}/\text{cm}^2$)
Scale Resolution	1ppm	1ppm	10ppm
Sampling Time	1.5 minutes	1.5 minutes	1.5 – 4 minutes
Tests per Box	4	5	5
Colour Change	Pink to White	Pink to White	Pink to White
Storage Conditions	25°C (77°F)	25°C (77°F)	25°C (77°F)
Kit Weight	367g (13oz)	250g (9oz)	208g (7oz)
Kit Dimensions	185 x 125 x 110mm	185 x 125 x 110mm	185 x 125 x 110mm
Part Numbers	E134----2	E134----1	E134----3

The Elcometer 134 units do not require the use of needles or contain mercury.

Elcometer 134A Chloride Ion Test Kit for Abrasives

Elcometer 134S Salt Detection Kit for Blast Cleaned Surfaces

Elcometer 134W Chloride Ion Test Kit for Water/Liquids

Test Method

The Elcometer 134 is simple and convenient to use with its three-step process

1. Pour CHLOR*EXTRACT™ solution into the latex tube.
2. Peel the protective backing off the flange of the tube, pinch to tube to retain the CHLOR*EXTRACT™ and stick the flange to the surface to be tested. Work the solution against the surface to extract the salts.
3. Peel the flange off. Insert the glass titration tube into the solution in the tube, and read the result.

