



elcometer® 2215 Lory LCH Cup



Conventional cylindrical cup with a needle fixed into the bottom for quick measurements on-site or during production.

The cup is first dipped into the product to be measured, then emptied. The flow time elapses when the point of the needle appears.

Part Numbers

Model	Cup Number	Range Centistokes (cSt)*	Part Number
Elcometer 2215/1	1	50-1100	K0002215M001
Accessories	Stopwatches		See Digital Stopwatch

Elcometer provide a viscosity recalibration service using state of the art, climate controlled facilities, for information please contact Elcometer

**for information only*



BAMR (Pty) Ltd, PO Box 23973, Claremont, 7735, South Africa
 Ph : 27 (0)21 683 2100, Fax : 27 (0)21 674 1485
 Email : sales@bamr.co.za, Web : www.bamr.co.za